


# U.S. Pharmacopeia method

## Abacavir Sulfate :Enantiomeric Purity



Column : CHIRALPAK® AD 0.46cmΦ × 25cmL (L51)  
 Mobile phase : See the gradient table  
 Solution A : Heptane / 2-Propanol / Diethylamine = 850 / 150 / 1 (v / v / v)  
 Solution B : Heptane / 2-Propanol = 1 / 1 (v / v)  
 Injection volume : 20μL  
 Column temperature : 30°C  
 UV detection : 286nm  
 (gradient table)


Time(min)	Solution A (%)	Solution B (%)	Flow Rate (mL/min.)
0	100	0	1.0
25	100	0	1.0
27	0	100	0.8
37	0	100	0.8
39	100	0	1.0
55	100	0	1.0

### System suitability

Sample: Transfer a quantity of USP Abacavir Stereoisomers Mixture RS to a suitable volumetric flask, add a volume of *Diluent* equivalent to 30% of the final volume, and sonicate until the solid is fully dissolved. Add a volume of 2-propanol equivalent to about 30% of the final volume, mix, and dilute with heptane to volume to obtain 0.4 mg/mL of USP Abacavir Stereoisomers Mixture RS.

Dilute: Methanol and trifluoroacetic acid (200:1)

Relative retention times: *trans*-Abacavir, Abacavir enantiomer, and Abacavir are 0.8, 0.9 and 1.0, respectively.



		Requirement	Result
Resolution	between <i>trans</i> -Abacavir and Abacavir enantiomer	≥ 1.0	1.7
	between Abacavir enantiomer and Abacavir	≥ 1.5	2.6