


## APPLICATION RESULTS


## COMPOUND IDENTIFICATION

Compound Name

Serine



## CHROMATOGRAM



## ANALYTICAL REPORT

Column Trade Name Column Size / Particle Size	CROWNPAK® CR-I(+) 150 x 3mm / 5µm
Mobile Phase Composition	perchloric acid aq. pH1.0 / acetonitrile = 85 / 15
Chromatographic Mode	Reverse Phase

Flow rate (ml/min)	0.1	Injection amount (mg)	0.002
Temperature (°C)	25	Pressure (MPa)	
Detection	200	Elution order	

No.	Rt	K'	Alpha Value	Resolution
1			1 - 2	