



APPLICATION RESULTS

COMPOUND IDENTIFICATION

Compound Name	Sertraline tetralone
---------------	----------------------



CHROMATOGRAM



ANALYTICAL REPORT

Column Trade Name Column Size / Particle Size	CHIRALPAK® AZ-H 250 x 4.6mm / 5µm
Mobile Phase Composition	n-hexane / ethanol = 80 / 20
Chromatographic Mode	Normal Phase Coated CSP

Flow rate (ml/min)	1	Injection amount (mg)	0.005
Temperature (°C)	25	Pressure (MPa)	
Detection	UV 254 nm	Elution order	

No.	Rt	K'	Alpha Value	Resolution
1	6.51	1.19	1 - 2 1.23	3.68
2	7.32	1.47		